
C-100 Circular Coil

- The coil is suitable for general-purpose stimulation.
- Equipped with power control and trigger button to support clinical operation.



Mechanical Properties

Weight of transducer head	0.6kg
Dimensions of transducer head	ø123 x 11.5 mm
Cable length	1.7m

Coil Winding Data

Inner diameter	20mm
Outer diameter	110mm
Winding height	6mm
Number of windings	14

Magnetic and Electrical Properties

Max initial dB/dt	35 kT/s near the coil surface.
Active pulse width	280µs (Biphasic)

Performance

Number of stimulations before warm-up at ambient temperature 20°C:	400 pulses
Mean output 75% of maximum at 1pps.	

Ordering Number	9016E058-
------------------------	-----------