
C-B60 Butterfly Coil

- The coil is suitable for focused stimulations.
- Equipped with power control and trigger button to support clinical operation.



Mechanical Properties

Weight of transducer head	0.7kg
Cable length	1.7m
Dimensions of transducer head	165 x 85 x 19 mm

Coil Winding Data

Inner diameter	35mm
Outer diameter	75mm
Winding height	11mm
Number of windings	2 x 10

Magnetic and Electrical Properties

Max initial dB/dt	35 kT/s near the coil surface
Active pulse width	280 μ s (Biphasic)

Performance

Number of stimulations, before warm-up at ambient temperature 20°C:	350 pulses
Mean output 75% of maximum at 1pps.	

Ordering Number	9016E048-
------------------------	-----------