



Medtronic

MAGPRO

60pps Option

Introduction

The MagPro 60pps Option is a software option for the MagPro R30 for the possibility for a higher stimulation repetition rate.

This document describes the technical data for 60pps Option which is different from the standard MagPro R30.

Please read the MagPro Family User Guide before this document.

Setting the Scale

MagPro R30 with 60pps Option

Repetition Rate settings

Timing	
Rep Rate	0.1, 0.2, 0.3, 0.4...1, 2, 3, 4, 5, 6...60pps
Trigger	
Rep Rate (ext. trig)	0.1, 0.2, 0.3, 0.4...1, 2, 3, 4, 5, 6...30pps



UK-edition



MagPro and Accessories are manufactured for/by:



DK Office

Medtronic A/S
Tonsbakken 16-18
DK-2740 Skovlunde
Denmark
Telephone: +45 44 57 90 00
Fax: +45 44 57 90 10
www.medtronic.com

Service Centers

Canada +1 905 826 6020 **Italy** +39 0266 1641 **Germany** +49 211 5293 211 **Austria** +43 1 240 440
Switzerland +41 21 803 80 00 **France** +33 1 5538 1700 **Belgium** +32 2 456 09 10 **Sweden** +46 8 760 83 00
ASEAN +65 677 0971 **Asian-Pacific** +852 2891 4456 **Middle East/Africa/South Asia** +961 1 370 670

Issued in Denmark, May 2007. P/N: M501-0291 (UK), rev. 1.0