

---

## MC-125 Circular Coil

- The coil is suitable coil for general purpose stimulation.
- Equipped with trigger button to support clinical operation.



### Mechanical Properties

Weight of transducer head	0.6kg
Cable length	1.3m
Dimensions of transducer head	ø130 x 11.3 mm

### Coil Winding Data

Inner diameter	28mm
Outer diameter	114mm
Winding height	6mm
Number of windings	13

### Magnetic and Electrical Properties

Max initial dB/dt	41 kT/s near the coil surface
Active pulse width	280μs (Biphasic)

### Performance

Number of stimulations, before warm-up at ambient temperature 20°C:	450 pulses
Mean output 75% of maximum at 1pps.	

<b>Ordering Number</b>	9016E055-
------------------------	-----------