
MC-P-B70 Placebo Butterfly Coil

- The coil's magnetic shield provides a field reduction of approximately 80%.
- The Placebo Coil has a mechanical outline and sound level identical to MC-B70.
- Number of stimulations before warm-up is identical to MC-B70
- Equipped with trigger button to support clinical operation.



Mechanical Properties

Weight of transducer head	1.8kg
Cable length	1.3m
Dimensions of transducer head	169 x 112 x 53 mm

Coil Winding Data

Inner diameter	20mm
Outer diameter	100mm
Winding height	6mm
Number of windings	2 x 10

Performance

Number of stimulations before warm-up at ambient temperature 20°C:	400 pulses
Mean output 75% of maximum at 1pps.	

Ordering Number	9016E059-
------------------------	-----------