

---

## MCF-P-B65 Placebo Butterfly Coil

- The coil's magnetic shield provides a field reduction of approximately 80%.
- The Placebo Coil has a mechanical outline and sound level identical to MCF-B65.
- With a reduction of stimulus intensity with 20-25% the coil can perform same number of stimulations as with the normal MCF-B65.
- Equipped with trigger button to support clinical operation.



### Mechanical Properties

Weight of transducer head	2.9kg
Cable length	2m
Dimensions of transducer head	174 x 94 x 53 mm

### Coil Winding Data

Inner diameter	35mm
Outer diameter	75mm
Winding height	12mm
Number of windings	2x (2x5)

### Performance

Number of stimulations before warm-up at ambient temperature 20°C:	1600 pulses
Mean output 75% of maximum at 1pps.	

<b>Ordering Number</b>	9016E060-
------------------------	-----------