
MagPro Remote Control



- The MagPro Remote Control provides the user with the possibility to operate the MagPro from a distance.
- It is by use of the MagPro Remote Control possible to make single stimuli and to set the output power from a distance.
- The LED indicates whether the Stimulator is enabled or disabled.
- Can be used on all available types of MagPro.

Mechanical Properties

Dimensions (WxDxH)	117 x 79 x 19 mm
Weight	0,4 kg
Cable length	8 m *
Connector	6 pole Lemo type
Encapsulation material	ABS plastic with soft side grip for hand-held comfort
Housing protection	IP20
Ordering Number	101-2211

* Other lengths can be specified.