

RT-120 Racetrack Coil



- The coil is elliptic in shape and is especially suitable for stimulation of wider areas such as bigger muscles.
- Equipped with trigger button to support clinical operation.

Mechanical Properties

Weight of transducer head	1.3 kg
Cable length	1.5 m
Dimensions of transducer head	ø90 x 200 x 26 mm

Coil Winding Data

Outer loop	ø80 x 160 mm
Inner loop	ø30 x 110 mm
Winding height	15 mm
Number of windings	10

Magnetic and Electrical Properties

Max initial dB/dt	31 kT/s near the coil surface
Active pulse width	280µs (Biphasic)

Performance

Number of stimulations before warm-up at ambient temperature 20°C:	1500 pulses
Mean output 75% of maximum at 1pps.	

Ordering Number	9016E064-
-----------------	-----------