

Sham Noise Generator



Technical Data

Output sham noise amplitude	Max 100dB
Pulse width of the sham noise	25-200mS

Mechanical Properties

Dimensions (WxDxH)	86 x 35 x 170 mm
Weight	365g
Cable length	2.5m
Connectors	Stereo MP3 input max 1V-rms 3.5 mm stereo jack Two Stereo audio outputs 3.5 mm stereo jack
Encapsulation material	ABS plastic
Housing protection	IP20 Tight

Ordering Number

Sham Noise Generator	9016C077-
Additional Headset	9016C078-

- In order to hide the click noise when a magnetic stimulation pulse is fired, a pulse of sham noise is sent into the ears of the patient.
- This sham noise pulse will hide the click noise from the coil for the patient; even at 100% stimulus intensity.
- For double blinded research experiments the MagPro operator also should receive the sham noise pulse.
- One headset included