



D-B80

Butterfly (figure 8), No cooling, No Placebo

Opened butterfly design for powerful stimulation. The coil is suitable for deep stimulation. The coil is slightly bent to closely follow curved shapes. Equipped with trigger button to support easy operation. Usable with a MagPro Compact stimulator with converter.

Contact

Website: www.magventure.com

Email: info@magventure.com

Performance

Number of stimulations before warm-up at ambient temperature 20°C: 500 pulses
Mean output 75% of maximum at 1pps.

Mechanical Properties

Weight of transducer head 0.9 kg/2 lbs
Cable length 1.7 m/5.58 ft.
Angle 120°

Magnetic and electrical properties

Max initial dB/dt 31 kT/s near the coil surface
Active pulse width 280µs (Biphasic)

Coil winding data

Inner diameter 67mm /2.64 in.
Outer diameter 95 mm /3.74 in.
Winding height 12 mm/0.47 in.
Number of windings 2 x (3+4)

Operating period

5 years

Ordering number

9016E043-