



MCF-P-B65

Butterfly (figure 8), Static cooling, Placebo

Placebo coil with mechanical outline and sound level identical to MCF-B65.
Equipped with trigger button to support easy operation.

Contact

Website: www.magventure.com
Email: info@magventure.com

Performance

Number of stimulations, before warm-up at ambient temperature 20°C (68°F): 1,600 pulses
Mean output 75% of maximum at 1pps.
Protocol: 60 trains @ 50 pulses/train @ 10pps @ Inter Train Interval: 25s @ Output=75% @ total number of stimulations 3,000

Mechanical Properties

Weight of transducer head	2.9 kg/6.4 lbs
Cable length	2 m/6.4 lbs
Dimensions of transducer head	174x94x53 mm /6.8x 3.7x2.1 in

Users in the USA

Caution: investigational device. Limited by Federal law to investigational use. In accordance with US Federal regulations, and IDE and/or IRB approval may be required.

Operating period

5 years

Ordering number

9016E060-