



MC-P-B70

Butterfly (figure 8), No cooling, Placebo

This is a placebo coil. The coil's magnetic shield provides a field reduction of approximately 80%

The coil has a mechanical outline and sound level similar to MC-B70. It has a slightly bend surface to closely follow curved shapes. Number of stimulations before warm-up is identical to MC-B70. The coil is equipped with a trigger button to support easy operation. Usable with MagPro Compact with converter.

Contact

Website: www.magventure.com

Email: info@magventure.com

Performance

Number of stimulations before warm-up at 400 pulses

ambient temperature 20°C/68°F

Mean output 75% of maximum at 1pps.

Mechanical Properties

Weight of transducer head 1.8 kg/ 4 lbs

Cable length 1.3 m/4.3 ft.

Dimensions of transducer head 169 x 112 x 36/53 mm (6.7 x 4.4 x 1.4/2.1 in.)

Angle 150°

Users in the USA

Caution: investigational device. Limited by Federal law to investigational use. In accordance with US Federal regulations, and IDE and/or IRB approval may be required.

Operating period

5 years

Ordering number

9016E0592