



TREATMENT CHAIR

The treatment chair may be used with the MagVenture TMS Therapy system. The special neck rest is designed for best possible comfort for the patient, offering a wide range of adjustment options. The neck rest in combination with a vacuum pillow will ensure stable positioning of the patient's head during treatment. The chair has 4 motors for height, tilting of seat, foot rest and back rest adjustment. All motors are controlled from a remote control. The chair comes with a foot rest cover and a terry chair cover.

Contact

Website: www.magventure.com
Email: info@magventure.com

Mechanical Data

Color	Grey (anthracite) PVC upholstery.
Width	63 cm/24.8 in without armrest, 80 cm/31.5 in with armrest
Height	63-87 cm (24.8-34.3 in)
Length	190-210 cm (74.8-82.7 in)
Weight	85 kg / 187 lbs
Patient max. weight	150 kg / 330 lbs
Safe working load	235 kg / 517 lbs

Electrical data

Mains inlet	230V~, 50/60Hz
Power consumption	Max. 50VA
Motors	4 motors for height, tilting of seat, footrest and backrest adjustment
Operating mode	Short-time duty cycle ED 10%: 1min ON, 9min OFF

Operating period

5 years

Ordering Numbers

9016B008-

9016B045- (Additional foot rest cover)

9016B046 -(Additional terry cloth cover)