



MARKING PLATES FOR DEPRESSION TREATMENT

With a marking plate mounted on the C-B60, the Cool-B70, the MCF-B65 or the MCF-B70, the DLPFC spot is easy to locate during the motor threshold determination. When APB motor cortex is located, a curved line is drawn with a pen along the marking plate. The curved line from the marking has the shape of the treatment coil which makes it easy to position the treatment coil on the scalp over the DLPFC.

Contact

Website: www.magventure.com
Email: info@magventure.com

Mechanical Data

Marking plate for C-B60 Coil	0.06 in. plastic. Based on 5.5 cm rule standard anterior to APB motor cortex
Marking plate for Cool-B70 coil	01-0.2 in plastic. Based on 5.5 cm rule standard anterior to APB motor cortex.
Clips-on marking plate for MCF-B65	0.06 in. plastic. Based on 5.5 cm rule standard anterior to APB motor cortex
Clip-on marking plate for MCF-B70	0.1 in. plastic. Based on 5.5 cm rule standard anterior to APB motor cortex

Ordering Numbers

- 9016B064- (Cool-B70 coil)
- 9016B050- (C-B60 coil)
- 9016B057- (MCF-B65)
- 9016B064- (MCF-B70)