



CAPS AND RULER

Treatment caps and measurement pin for use with the MagVenture TMS Therapy system when determining motor threshold. With the caps, it is easy to mark the placement position for the magnetic coil. This facilitates correct repositioning of the coil at future treatment sessions. Specific data pertaining to the individual patient can be written directly on the cap, such as patient ID and distance to the edge of the cap. To measure the distance, we offer a small measurement pin with scale in millimeters and centimeters. Optional chin straps are also available.

Contact

Website: www.magventure.com
Email: info@magventure.com

Mechanical Data

Caps	Material	Textile 97% cotton / 3% lycra with elastic band in the back of the neck
	Sizes	Extra small: 52-54 cm/20.5-21.3 in Small: 54-56 cm/21.3- 22 in Medium: 56-58 cm/22-22.9 in Large: 58-60 cm/22.9-23.6 in Extra large: 60-62 cm/23.6-24.4 in

Ordering Numbers

- 9016B059- (Textile Cap, Extra small – 10 pcs) Chin straps not included.
- 9016B020- (Textile Cap, Small – 10 pcs). Chin straps not included.
- 9016B021- (Textile Cap, Medium – 10 pcs). Chin straps not included.
- 9016B022- (Textile Cap, Large – 10 pcs.) Chin straps not included.
- 9016B023- (Textile Cap, Extra Large – 10 pcs.) Chin straps not included.
- 9016B060- (Textile Caps, bundle – 10 pcs. of XS-S-M-L-XL) Chin straps included.
- 9016B065- (Chin Straps for textile caps)– 10 pcs.)