



MAGVENTURE 360

The MagVenture 360 patient management system is a cloud-based, GDPR compliant system for use in a clinical setting. It is an optional add-on to the MagVenture TMS Therapy system. With this system you get anytime, anywhere access from all platforms such as laptops and tablets, safe data storage and the possibility to store all patient data in one place and cut down administration time. The solution consists of a 360 gateway unit with native application and a monitor/touch screen with data from the 360 web application along with a computer (not supplied by MagVenture) with internet connection.

Contact

Website: www.magventure.com
Email: info@magventure.com

Mechanical Data

360 Gateway

Dimensions (HxWxD): 184 x 140 x 40mm/7.2 x 5.5 x 1.6 in
Weight 0.95 kg/2.1 lbs

360 Monitor

Dimensions (HxWxD): 250x335x29 mm/9.8 x 13.2 x 1.2 in
Weight 2.5 kg/5.5 lbs
Display and touch screen 12.1 inch. 1024x768

Environmental Data

Operating Temperature: 10-30°C/50-86°F
Storage Temperature: 5-50°/41-122°F
Operating Humidity: 30-60% RH
Storage Humidity: 20-80% RH
Operating Pressure: 80-106kPa
Operating Altitude: -400 m to 2.000 m/-1312ft. to 6561ft.
Storage/transport Pressure: 24-106Pa
Storage/transport Altitude: -400m to 10.000 m/-1312ft. to 32808ft.

Operating period

8 years

Ordering Numbers

9016S019- (360 software – annual license)
9016B070-(360 Mounting Kit for Monitor)
9016C090- (360 Gateway unit)
9016C091- (360 Monitor)