



COIL COOLER UNIT

The coil cooler unit is used as external cooling for cool coils. It is equipped with special liquid cooling media. Can be used in combination with the high-performance option for coil cooler to allow for even more aggressive protocols.

Contact

Website: www.magventure.com

Email: info@magventure.com

Mechanical Data

Weight	10 kg/22 lbs
Dimensions (HxWxD)	200x300x300 mm (7.9 x 11.8 x 11.8 in)
Mains Voltage:	230V~, 50/60 Hz
Power consumptions:	Maximum: 40VA
Fuses F1, F2:	T1.6A/250V~, IEC 60127-2/3

Ordering Numbers

9016B015-