



TROLLEY FOR MAGPRO X/R-MODELS

The trolley is designed to hold a complete system with stimulator, isolation transformer, coil cooler unit, and vacuum pump unit. Prepared for mounting of accessories for trolley, flexible arm and sham noise generator. The trolley comes in both a short and a high model – the high trolley has an extra shelf on the top for navigation systems or customer's external systems, for example. Both the high and the short trolley is suitable for MagPro R30, MagPro R30 with MagOption, MagPro X100 and MagPro X100 with MagOption and MagPro R100.

Contact

Website: www.magventure.com

Email: info@magventure.com

Mechanical Data

Short version:

Weight:

17kg/37.5 lbs

Dimensions (HxWxD):

80x64x55 cm/ 31.5x 25.2x21.65 in

High version:

Weight:

25kg/55 lbs

Dimensions (HxWxD):

128x64x55cm/50.39x25.2x21x65 in

Ordering Numbers

9016B010- (short trolley)

9016B043-(high trolley)