



## COIL HOLDERS FOR TROLLEYS

MagVenture offers coil holders for all our different types of coils whether it is coils without cooling, coils with static cooling or with active cooling. The coil holders can be attached on either the trolley or on a wall.

### Contact

Website: [www.magventure.com](http://www.magventure.com)

Email: [info@magventure.com](mailto:info@magventure.com)

### Technical data

Universal holder for coils $\varnothing 25$ - $\varnothing 38$ :	Fits all coils
Holder for coils with $\varnothing 25/\varnothing 29$ handle:	For coils without cooling and older static cooled coils
Holder for coils with $\varnothing 29/\varnothing 38$ :	For newer static cooled coils and active cooled coils

### Operating period

N/A

### Ordering Numbers

9016B036- (For all coils)

9016B032- (Holder for coils with  $\varnothing 25/\varnothing 29$  handle)

9016B035- (Holder for coils with  $\varnothing 29/\varnothing 38$  handle)