



## **ELECTRODE CABLE**

Shielded electrode cable for MEP Monitor for active, reference and ground electrodes. Available with two different connector types.

### Contact

Website: [www.magventure.com](http://www.magventure.com)

Email: [info@magventure.com](mailto:info@magventure.com)

### Technical data

#### Shielded Electrode cable

Cable length	3 m/9.8 ft
Connector for MEP Monitor	5-pole DIN 240°
Connector for electrodes	1.5mm/0.06 in male TP (3 pcs.)
Ordering Number	9016C081-

#### Shielded Electrode cable

Cable length	2 m/6.6 ft
Connector for MEP Monitor	5-pole DIN 240°
Connector for electrodes	0.7mm/0.03 male TP (3 pcs.)
Ordering Number	9016C085-