



EXTERNAL CONTROL SOFTWARE FOR MAGPRO

If the MagPro stimulator is placed in an equipment room behind the MRi-scanner room and away from the operator room, this option can be used to externally control the MagPro and receive log information. With the external control software, it is possible to control different settings on MagPro from a standard PC with USB or Serial port interface.

Contact

Website: www.magventure.com

Email: info@magventure.com

Technical data

Stimulator compatibility: MagPro R100, MagPro R30, MagPro R30 with MagOption, MagPro X100 and MagPro X100 with MagOption.

Operating period

N/A

Ordering Numbers

9016S014-